



Report No:ATCCR20113011F

Test Report

Sample Category Spray Dressing

Client Beijing Gold Mountainriver Tech Development Co.,Ltd.

Test Category Test Entrust

Date of Report 2020.12.12





Detection Information

Client	Beijing Gold Mountainriver Tech Development Co.,Ltd.		Sample Source	Inspect
Client address	No.30 Baishiqiao Road, Haidian District, Beijing 100081,P.R. China		Sample State	Normal
Date of Receives samples	2020.11.30		Date(s) of tests	2020.11.30-2020.12.11
Sample No	ATCCR20113011F-1130CP01			
Category	Test Project	Test Standard and Method		Test Instruments
Spray Dressing	Virus inactivation SARS-CoV-2	ISO 18184:2019 Textiles Determination of antiviral		Biosafety Cabinet

End

Remarks: Beijing Gold Mountainriver Tech Development Co. Ltd

Report Preparer: 
 Report Reviewer: 


Authorized Signatory: 
 Date of Issues Report: 2020.12.11
 (Special Chapter for Inspection and Inspection)



Test results

Virus Types	NO	$\lg(V_{a0h})$ ($\lg\text{TCID}_{50}/\text{mL}$)	$\lg(V_{b4h})$ ($\lg\text{TCID}_{50}/\text{mL}$)	$\lg(V_{c4h})$ ($\lg\text{TCID}_{50}/\text{mL}$)
SARS-CoV-2 MDCK cells	1	6.68	6.57	4.25
	2	6.66	6.52	4.23
	3	6.69	6.54	4.21
Average Value of $\lg\text{TCID}_{50}/\text{mL}$		6.68	6.54	4.23
Antiviral Activity Value		2.31		
Antiviral Activity Rate (%)		99.51		

Sample Picture:



End

天
天
天